


StethViewtm

Digital Recording Software

Auscultation - the key to physical diagnosis

at the point of care




**StethView Two-Channel PC Display
Heart Sounds and ECG**



- 📌 **Display Heart and Breath Sounds** - View sounds on a PC as a phonocardiogram or phonopneumogram.
- 📌 **Record and Play Back Sounds** - Play back recorded sounds at normal or half-speed.
- 📌 **One or Two Channel Display** - Permits simultaneous viewing of phonocardiogram and ECG or phonocardiogram and phonopneumogram.
- 📌 **Installed Database of Simulated Sounds** - Permits teaching at the bedside with a variety of sounds
- 📌 **Compare Sounds** - Compare one sound with another for teaching purposes.
- 📌 **General Clinical and Teaching Use** - Used with Cardionics Auscultation Systems listed on reverse side.



StethView Display



With **StethView** Software installed on a current PC, the patient is auscultated as usual using any of these Cardionics Auscultation Systems:

- The E-Scope Electronic Stethoscope
- The Ausculette Bedside Auscultation System
- The SimulScope Bedside Auscultation System with the Wireless HeartMan PC Receiver
- The Classroom Infrared Emitter
- The Auditorium Infrared Sound System

Sounds are seen on the PC monitor as a phonocardiogram (PCG). Sounds can be recorded up to 120 seconds and playback quickly and easily. The user has a choice of a one or two channel display. Two channels can be used for heart sounds and ECG, if desired. Sounds can be edited as desired or exported as a .wav file.

StethView has a databank of simulated sounds for teaching purposes. The sounds supplied are: S3, S4, MVP, Aortic Stenosis, Aortic Regurgitation, Mitral Stenosis, Mitral Regurgitation, Pericardial Rub, Hypertrophic Cardiomyopathy, Pulmonary Hypertension, Summation Sound, VSD, PDA, ASD, and Innocent Systolic Murmurs.

With Cardionics Auscultation Systems and the computing power of **StethView** Software, we're revolutionizing auscultation - right before your ears! Give us a call - We'll arrange a demonstration!

Cardionics, Inc.

910 Bay Star Blvd. Webster, Texas 77598 USA
Telephone : 281-488-5901 US and Canada : 1-800-364-5901
Fax: 281-488-3195 Email: keithj@cardionics.com Web site: www.cardionics.com